HW-503-A

A multi-function surgical operation table, which is powered by electric mechanical transmission and safe voltage control. Surgical operation table is mainly used for orthopaedics, urology, anorectal, chest, gynecology, and obstetrics operation. Sliding movement is suitable for C-arm operation.



ltem No.	HW-503-A
Dimension	2000x500mm
Elevation	750-1050mm
Back Section	Up 65°/Down15°
Kidney Bridge	120mm
Trendelenburg	+25° / - 15°
Lateral Tilt	Left 18°/Right18°
Sliding Movement	300mm
Leg Section	Down 90° / Outreach90°

