

## Electric home care Bed

Bed size	Load capacity	High-low
2080x960x500mm	≤ 250kg	550-740mm
Back raised	Leg raised	Left-turn&right turning
0-75°	30-45°	0-45°



PE-E02

### Specification:

- 1)1.5mm bed frame and 1.2mm bed platform
- 2)ABS head and foot board
- 3)Aluminium guard rails
- 4)Built-in commode:215\*250mm

Backrest up range:0-75° (+5° )

Leg down & up: 30-45° (+5° )

Left-turn&right turning : 0-45° (+5° )

Safe load:250kg

Overall Approx Dim(mm):2080x960x500mm

