

WYK12 Ceiling Shadow less Operating Lamps



WYK series shadowless operation lamps with good performance in illumination, flexible and stable suspension systems contents varies requirement on height and angle. They are suitable for the requirements of various operation conditions. It can be used as sole illumination for minor operation and auxiliary illumination for major operation.

- Excellent bright cool light
- Smooth adjustment of light intensity
- Smooth regulation of light beams focusing
- With the ultra-thin optical lens
- Illumination depth up to 600mm
- Active system of shadow reduction
- Perfect fidelity of color rendering of illuminated objects
- Mobility and an excellent movement range of the light head

Technical Parameter	Model No.
	WYK12
Intensity (Lx)	≥100,000
Color Temperature (K)	3800±500
Color rendering index (Ra%)	≥80
Light Field Diameter (mm)	150±35
Depth of Illuminance (mm)	600
Input power (Watts)	150va
Bulb power (Watts)	300
Brightness	Fixed/ Adjustable
Diameter of lamp head (mm)	900