

CE

WYK9

Ceiling Shadowless Operating Light



User-friendly and Best Value for Money



CAD/CAM designed multi-angle reflectors, eliminate shadows from the surgeons heads and hands



Switch on/off



Stable suspension system



Light focusing adjustable



Brightness adjustable

WYK9

Ceiling Shadowless Operating Light

Premium Quality, Surprising Price

WYK9

Ceiling Shadowless Operating Light

General Features

Good performance in illumination depth, CRI;

Flexible and stable suspension system;

Excellent shadow control;

Best value for money.

Technical Specifications

Illumination: 40,000-90,000 lux

Color temperature: 3800±500K

Color rendering index: 80%

Light field size: 150±35mm


Illumination depth: 450mm

Bulbd power: 24V/25W

Power: 225W

Diameter of lamp head: 900mm

Temperature rise at surgeon's head: Less than 5°C



All measurement is in millimeter.