


上海卫远医疗器械有限公司 SHANGHAI WEYUAN MEDICAL DEVICE CO., LTD



EXCEED IN QUALITY AND MORE

WYLEDK4

LED Operating Light





上海卫远医疗器械有限公司 SHANGHAI WEYUAN MEDICAL DEVICE CO., LTD



EXCEED IN QUALITY AND MORE

WYLEDK series surgical lighting system is designed to illuminate operating room surgical field and can be adjusted for light intensity in a variable pattern size.

General Features

High-intensity of 80,000 lux at 1 meter;

4 500K color temperature offers natural and true coloring;

Accurate color rendition;

Super deep-cavity illumination;

Excellent shadow control;

Energy efficiency;

Osram LED chips – each rated with a 50,000 hour bulb life;

Extremely smooth focus mechanism with the focusing knob;

Variety of configuration: floor standing, wall mounted and ceiling.

Technical Details

Illumination: 80.000 lux


Color temperature: 4000±500K Color rendering index: 90% Light field size: 200±50mm Illumination depth: 10000mm

Bulb power:1-3W

Diameter of lamp head: 500mm

Qty of LEDs: 20

Temperature rise at surgeon's head: Less than 1°C



Measurement: mm