

CE

EXCEED IN QUALITY AND MORE



# **LED Operating Light**





## 上海卫远医疗器械有限公司 SHANGHAI WEYUAN MEDICAL DEVICE CO., LTD



#### EXCEED IN QUALITY AND MORE

WYLEDK series surgical lighting system is designed to illuminate operating room surgical field and can be adjusted for light intensity in a variable pattern size.

## **General Features**

High-intensity of 130,000 lux at 1 meter;

4 500K color temperature offers natural and true coloring;

Accurate color rendition;

Super deep-cavity illumination;

Excellent shadow control;


Energy efficiency;

Osram LED chips – each rated with a 50,000 hour bulb life;

Extremely smooth focus mechanism with the focusing knob; Variety of configuration: floor standing, wall mounted and ceiling.

## **Technical Details**

Illumination: 130,000 lux Color temperature: 4000±500K Color rendering index: 96% Light field size: 200±50mm Illumination depth: 10000mm Bulb power:1W Diameter of lamp head: 500mm Qty of LEDs: 25 Temperature rise at surgeon's head: Less than 2°C



Measurement: mm