

CE

WYLEDK9

Ceiling LED Operating Light



User-friendly and Best Value for Money



Osram LED lighting deliver cool, white bright light;


CAD/CAM designed multi-angle reflectors, eliminate shadows from the surgeons heads and hands



Switch on/off & brightness adjustment



Stable suspension system



Light focusing adjustable

WYLEDK9

Ceiling LED Operating Light

Premium Quality, Surprising Price

WYLEDK9


Ceiling LED Operating Light

General Features

- High intensity
- Cool, bright white light
- Accurate color rendition
- Superior deep-cavity illumination
- Excellent shadow control
- Energy efficiency

Technical Specifications

- Illumination: 40,000-150,000 lux
- Color temperature: 4000±500K
- Color rendering index: 90%
- Light field size: 200±50mm
- Illumination depth: 1000mm
- Qty of LEDs: 45pcs
- Power: 1W*45
- Diameter of lamp head: 1000mm
- Temperature rise at surgeon's head: Less than 1°C



All measurement is in millimeter.