## Vital signs monitor



### **Features**

- Compact and flexible appearance
- Simple and convenient operation, which is applicable for community medical treatment and family ect.
- 4.3 inch color TFT display and highlight LED display Adjustable audible and visual alarms
- Built-in rechargeable lithium battery for uninterrupted
- monitoring

### **NIBP**

- Method: Oscillometry
- Measure mode: Manual, Auto, STAT
- Unit: mmHg, kPa
- Measure and alarm range:
- --- Adult Mode

SYS 40 ~ 280 mmHg DIA 10 ~ 225 mmHg MEAN 20 ~ 240 mmHg

Pediatric Mode

SYS 40 ~ 220 mmHg DIA 10 ~ 160 mmHg MEAN 20 ~ 170 mmHg

Neonatal Mode

SYS 40 ~ 135 mmHg DIA 10 ~ 100 mmHg MEAN 20 ~ 110 mmHg

- Resolution: 1mmHg
- Accuracy: ±5mmHg

# **Technique Specification**

#### SPO2

• Measuring Range: 0 ~ 100 %

Resolution: 1%

Accuracy: 70% ~ 100% (±2 %)

Pulse Rate: 20-250 BPM

Resolution: 1 BPM

Accuracy: ±3 BPM

## TEMP

- Measure and Alarm Range: 0 ~ 50 °C
- ◆ Resolution: 0.1°C
- Accuracy: ±0.1°C

### **Standard Parameters:**

NIBP, SPO2, PR, TEMP

## Packing:

- → Single package size: 35 \*25 \*26cm
- Single gross weight: 3 KG
- 8 unit per carton, package size: 72.5 \*53
  \*56.5 cm, total gross weight 24 KG