



ABS刹车系统  
ABS braking system



智能方向操作杆  
Intelligent manual  
steering lever



手电动模式切换  
Hand power mode  
switch



车架快速折叠  
Frame lock quick



加大防滑轮胎  
Increase the  
anti-slip tyre



MODEL :

**A2 PLUS RECLINING**

车长	Overall length	110cm	载重量	load weight	120公斤
车宽	Vehicle Width	63cm	爬坡能力	Climbing Ability	≤ 13°
车高	Overall height	120cm(含枕头)	电机功率	The Motor Power	24V DC250W*2
座宽	Base width	45cm	蓄电池	Battery	24V12AH/24V20AH
前/后轮尺寸	The Front / Rear Wheel Size	8/12寸	续航里程	Range	10-20KM(视路况和电池容量)
整车净重	The Vehicle Weight	29kg+ 3kg(锂电池)	时速	Per Hour	1-6KM(五档可调)