



EABS刹车系统
EABS braking system



智能方向操作杆
intelligent directional
controlling lever



手动模式切换
Hand control mode
switch



车架快速折叠
Frame fast folding



加大防滑轮胎
Increase the
anti-slip tyre



后减震



侧面调节靠背



MODEL :

DELUXE SPORTS

车长	Overall length	122cm	载重量	load weight	120公斤
车宽	Vehicle Width	65cm	爬坡能力	Climbing Ability	≤ 13°
车高	Overall height	128cm	电机功率	The Motor Power	24V DC250W*2
座宽	Base width	45cm	蓄电池	Battery	24V12AH/24V20AH
前/后轮尺寸	The Front / Rear Wheel Size	10/16寸	续航里程	Range	10- 20KM(视路况和电池容量)
整车净重	The Vehicle Weight	41kg+ 10kg(电池)	时速	Per Hour	1-6KM(五档可调)