



thin & light Notebook design

SS-350

Full Digital Ultrasonic
Diagnostic System



Notebook *Arise*

SS-350

Full Digital Ultrasonic Diagnostic System



Feature:

- 12-inch LED monitor;
- Abundant presets with one key optimization;
- Software packages: abdomen, gynecology, obstetrics, small parts, vessels, MSK, cardiology, etc.;
- Report template includes images, measurement data and summary;
- Friendly interface for easy and fast operating progress;
- Built in Image storage : ≥ 64 frame. External USB storage, less than 9 seconds (the storage speed is up to USB device.);
- Body marks: ≥ 97 kinds;
- With rechargeable lithium-ion battery, it can be standby for about 5 hours (optional);
- Biopsy guidance is also available.

Standard Configuration:

- Main Unit, 3.5Mhz Convex probe

Options:

- 6.5MHz transvaginal probe, 7. 5MHz linear probe, video printer, rechargeable lithium-ion battery

