



418-0012



418E0012



Fixation Cards

3.5V Streak Retinoscope

3.5V Streak Retinoscope is designed exclusively to determination of the refraction of the eye. The streak retinoscope is found by most practitioners to be easy to use, fast, accurate, and especially valuable in determining the axis of astigmatism.

- ◆ Multi-coated optics
- ◆ A single control sleeve for vergence and rotation. Comfortable operation.
- ◆ Totally dustproof, maintenance free
- ◆ The lowest position gives a divergent beam.
- ◆ Superior streak quality enables easier and quicker refractions—even in patients with small pupils
- ◆ The powerful, bright bulb ensures easy recognition of neutralisation.
- ◆ Fixation Cards for dynamic retinoscopy

418-0012 3.5V Streak Retinoscope with a halogen bulb, fixation cards, without handle

418E0012 3.5V LED Streak Retinoscope with LED bulb, fixation cards, without handle

Accessories

418-0050 Fixation cards (set of 4)

418-0910 R4180 3.5V Halogen bulb

418E0910 R4180-LED 3.5V LED bulb

Order Code	3.5V Streak Retinoscope	3.5V Streak Retinoscope	3.5V Handles				Power Source			Case	
	Halogen	LED	U-100	U-100	U-100	C-size	DT801 / DT803	486-1000	486-0010	Carton box	Zipper case
418-1811	■		■				★				■
418-1111	■		■				★			■	
418-5811	■			■			★			■	■
418-5111	■			■			★			■	
418-4812	■				■						■
418-2002	■							■		■	
418-3002	■								■	■	
418E1811		■	■				★				■
418E1111		■	■				★			■	
418E5811		■		■			★			■	■
418E5111		■		■			★			■	
418E4812		■			■						■
418E6812		■				■					■
418E2002		■						■		■	
418E3002		■							■	■	

★ The kit does not include the charger. Pls refer page 26 power sources to order DT series charger separately.