



Nasal Mask



Air filter
× 5



User manual



Air tubing



T-junction



Power supply



Travel bag


Auto BIPAP VM-8

American Sensor

Intelligent Breathing Algorithm

Improve Overall Breathing

German Technology



Patient-device synchrony
Real-time displaying



Whisper quiet
light and compact



Efficient heating
humidifier system



Professional sleep study report

Unique sleep study report of whole week & month.

New operating system is easier to operate, interface is more exquisite.



Specification

Description	Auto BIPAP VM-8
Work Mode	CPAP, Auto,S, T, ST, Auto S, APCV
Sound	<28 dB
Screen	2.8 inch TFT, 320 x 240
Pressure Range	4-30 cmH2O
IPAP	4-30 cmH2O
EPAP	4-20 cmH2O
Power Supply	Input 100-240V, 50-60Hz, 1-2A; Output + 24V, 2.5A
DC Voltage	24VDC, 2.5A
Over-Temperature Protection	Stop automatically, alarm with 4 times
Lock Function	Yes
Screen Auto Sleep	Auto sleep in 15 minutes, press any button to highlight screen
Backup Rate	5-50min, adjustable
I/E	10-80%
Ramp Item	0-60min adjustable 1 min/Step
Auto ON/Off	Yes
Inspirate Triggler	1-5 level
Expirate Level	1-5 level
Inspirate Sensitivity	1-5 level
EPR	Yes, 0-3 level
Auto-Altitude	Yes
Leak Compensation	Yes
Awake Detection	Yes
Breath Event Detection	Yes
Water Capacity	200ML max
Humidifier Level	Integrated Heated, 0-5 level
Alert	Power off, mask off
Dimensions(cm)	29L x 15W x 10H(machine)
Carton Dimensions(cm)	34.5L x 20.5W x 25.5H(1pc carton)
	40L x 34.5W x 25.5H(2pcs carton)
	59.5L x 34.5W x 25.5H(3pcs carton)
Weight(with humidifier)	1.7Kg(net weight); 4.0Kg(gross weight)
Work For	OSA,CSA,MSA,COPD



Hunan VentMed Medical Technology Co., Ltd.

Factory: Baoqing industrial park, Shaoyang, Hunan, China

Office: +86-739-5188959 Web: www.ventmed.net

