



Handheld oximeter

P01



Features

- Compact and potable design
- Numeric display with plethysmogram display
- 1.77 inch color TFT LCD in real time display, displayable in big front and big screen
- Adjustable audio and visual alarm
- Built-in Li-ion Battery for up to 8 hours continuous working

Standard Configuration:

- Oximeter main unit: 1 PC
- Adult finger SpO2 sensor: 1 PC
- USB communication cable: 1 PC
- Instruction manual: 1 PC
- Gift box: 1 PC

Packing:

Single package size: 16.5\*12.2\*7.2cm  
 Single gross weight: 0.25KG  
 32 unit per carton, package size: 48\*32\*36cm, total gross weight 10.5 KG

Optional accessories:

- TEMP sensor
- Pediatric finger/Neonatal Wrap SpO2 sensor
- Nellcor SpO2 sensor

Specification

- Parameter: Spo2, Pulse Rate
- Spo2 Range:

Range: 0-100%  
 Resolution: 1 %  
 Accuracy: ±2% at 70-99%  
 0% - 69%: Unspecified

- Pulse Rate:

Range: 30bpm ~250bpm  
 Resolution: 1bpm  
 Accuracy: ±2% at 30-250bpm

Application

- Handheld Pulse Oximeter is suitable for being used in family, hospital, oxygen bar and community health center etc.

Measuring parameter:

- Spo2
- PR