

18 Body Data Body Fat Analyzer H01A

Maximum Weight Recommendation	200KG
Health Scale Metal Bonding Point	Single-point Type
Function	Weight Measuring, Fat Measuring, Water Content Measuring, Body Fat and Water Content Testing
Display Type	Digital Display
Scale Use	height weight scale
Product name	Height and Weight Scale
Capacity	200kg
Unit	Kg Lb St
Display	LCD(LED)Display
Single package size:	105X50X30 cm
Single gross weight:	7.300 kg



* ABS plastic platform		
* LED Display size: 105*69 mm		
* Max weight: 200 kg, Height: (95cm -195cm)		
* Division: 50g , Height division: 0.5cm	Dimensions of the scale:	44.8 x 34 x 207.5 cm
* Step on technology for weight function	Carton box volume (1ps/CTN):	134 x 39x 22 cm
* Low battery /overload indication	Net weight:	9.5 kg
* Power : USB Charge and Storage the power	Gross weight	11 kg



① HD Voice Broadcast

② Can store test data

③ HD LCD electronic screen

④ Supports 5V DC power supply

At the moment, drinking a glass of water is a perceptive experience



The ultrasonic probe can easily record various data of the body and measure height/weight/fat at the same time And other data, can store 30 groups of personal data, can be searched records

- ✔ Ultrasonic measurement
- ✔ Voice broadcast
- ✔ 18 items of body data
- ✔ Bluetooth connection
- ✔ Free SDK

18 Body data



45.05 kg Weight	153.0 cm Height	19.2 BMI	36 Metabolic Age	8 Visceral fat Index
18.7% BFR	47.3% Muscle Rate	62.2% Body water	49.8 kg Standard weight	4.8 kg Weight control
2.7 kg Bone mass	1626 kcal BMR	18.6% Protein Rate	36.6 kg Weight Without fat	21.3 kg Muscle Mass