



lovego

Go along with love

Portable Oxygen Concentrator

SG02L

(Spring Gale™)



50%
more O₂

lovegopoc.com



50% More Oxygen Supply

SG is a new portable oxygen concentrator (POC) with a 50% larger supply of oxygen in pulse flow. Compared to traditional POC's 210ml per flow setting, SG supplies 333ml for each setting.



Traditional POC	210 ml	210 ml	210 ml	210 ml	210 ml	1050 ml
	1	2	3	4	5	flow setting
New POC (SG)	333 ml	333 ml	333 ml	333 ml	333 ml	1680 ml
	1	2	3	4	5	flow setting



Even at the same flow setting of 5, SG provides 50% more oxygen gas than other devices.

With a sufficient supply of oxygen, patients are capable of daily life and outside activities with less effort. It solves the problem of insufficient oxygen when patients need more oxygen in outside or in the future – a common complaint from previous small POCs.



Use with Confidence

SG has a balance of large flow oxygen production, is lightweight, small in size, and has a longer battery time. Its 4.3" large screen with a touch function, mobile phone App remote control, and durable designs fully consider patients' daily lives and outside activities.



4.3" Touch Screen

Lovego has a 4.3" screen with a wonderful touch design. It's two or three times bigger than other devices. A larger, clearly visible screen would be more friendly to the clients' daily lives.

Drop-resistant

SG is constructed from top-quality magnesium alloy, a metal material used in aerospace engineering, designed to withstand falls, ensuring long-term reliability.





Lovego


Lovego **App**

– Manipulate Device on the Mobile Phone

It breaks the barriers of physical distance and offering true freedom. You can control and manipulate the machine through the mobile phone directly

Everything is visible

You just need to touch your phone to control your machine when on the go, sitting in wheelchairs, and lying in bed. You don't have to struggle to get the machine and then press the button.



Technical Specifications



SG Highlights

- 50% more oxygen supply
- Handy App to life
- Durable to use
- Lightweight
- Long battery and fast charging

SG device

SG02L

Weight:	5 lbs (2.28 kg)
Dimensions: (L X W X H)	7.32' ' x 3.54' ' x 8.07' ' (18.6 cm x 9 cm x 20.5 cm)
Housing construction	Aerospace materials – Magnesium Alloys Resistant to falls and impacts, Flame-retardant Resistant to extreme cold and high temperatures

Oxygen delivery:	Pulse Flow
Flow settings:	1, 2, 3, 4, 5 (333ml per setting, 50% more oxygen)
Oxygen concentration:	90% (+5%/-3%)
Altitude:	0 to 16,400 ft 0 to 5000 meters

Battery Life:	Setting	Standard battery	Extra battery
	1	4 hrs	8 hrs
	2	3 hrs	6 hrs
	3	2 hrs	4 hrs
	4	1.5 hrs	3 hrs
	5	1.2 hrs	2.4 hrs

Charging time: 1 hour both in the car and home power

Sound: 43 – 55 dBA

Product Codes

SG Machine	SG02HM
Standard Battery	SG02IB
Extra Battery	SG02IB
AC Adaptor	SG02AA
DC Charger	SG02DC
Shoulder Bag	SG02SB
Nasal Cannula (2.5 meters)	SG25NC
User Manual	SG02UM

**SG standard package includes the SG portable oxygen concentrator, SG shoulder bag, one SG standard battery, one AC Adaptor, one DC Charger, one nasal cannula, and one user manual.



Lovego, go along with love!